



Universidad
del País Vasco
Euskal Herriko
Unibertsitatea

Inline monitoring of particle size in emulsion polymerization processes by Photon Density Wave (PDW) Spectroscopy

¹U.O. Aspiazu, ¹M. Paulis, ¹J.R. Leiza




NANOPAT:

¹ POLYMAT, Kimika Aplikatua saila, Kimika Fakultatea, University of the Basque Country UPV/EHU, Joxe Mari Korta zentroa, 20018 Donostia-San Sebastián.

Table of contents




- 1. Introduction**
 - i. NanoPAT
 - ii. PDW background
- 2. Analysed process**
- 3. Results**
- 4. Conclusions**



1 i) Introduction.

NanoPAT European Project: Application of three new real-time analytical tools for particle size (d_p) and particle size distribution (PSD) in-line or on-line monitoring.

Present Work objective: Assessing Photon Density Wave (PDW) spectroscopy analysis method as in-line monitoring technique in emulsion polymerization processes.



CryoTEM image for final latex ($t=130\text{min}$). Scale: white line in the base of the figure=500nm

1 ii) PDW background





Phase shift (k_ϕ)
Intensity shift (k_I)



Absorption Coefficient (μ_a)
Reduced Scattering Coefficient (μ'_s)



Particle diameter (d_p)
(Particle size distribution. PSD)



1 ii) PDW background

Phase shift (k_ϕ)
Intensity shift (k_I)



Absorption Coefficient (μ_a)
Reduced Scattering Coefficient (μ'_s)




Particle diameter (d_p)
(Particle size distribution. PSD)

$$\mu_s' = \frac{3\phi n Q_s}{2d_p} \int_{-1}^1 f(\cos(\theta)) S(\theta, \phi) [1 - \cos(\theta)] d\cos(\theta)$$


*Hass, Bressel, Münzberg, Reich, Appl. Opt. 52, p. 1423-1431 (2013)

1 ii) PDW background


2 potential solutions



Which is the
correct one?



Comparison of data
obtained with
different λ




2. Analysed process

All acrylic latex synthesis

In-situ seed formation
(30 min)

Particle growth
(100 min)

Postpolymerization
(30 min)





Recipe:

- Monomer (w/w)
 - MMA 51 %
 - BA 47 %
 - MAA 2 %
- Initiator: APS
- Surfactant: SDS
- Water

erriko
tatea

3. Results. Inline analysis




Effect of laser wavelength



3. Results. Inline analysis

Good agreement between inline PDW and offline DLS results

For the range 50-120 nm and up to 35 % SC:
the smaller the wavelengths, the better the inline analysis.



Analysis of broader dp and SC range

Coagulation detection

Bimodal latex analysis







Universidad
del País Vasco
Euskal Herriko
Unibertsitatea


POLYMAT

3. Results. Particle size range analysis

SC40 %; monomodal; dp $175 \rightarrow 400$ nm





3. Results. Coagulation detection








Feeding rate: 1 g/min
Feeding time: 65 min
Final solution concentration: 0.05-0.15 M

* Wavelength used for the analysis: $\lambda = 636 \text{ nm}$



3. Results. Detection of bimodal latex formation

SC 40 %; bimodal; dp 45→80/300 nm (Post-process)



4. Conclusions

- ✓ Accurate monitoring of particle size during seeded semibatch emulsion polymerization processes (SC of 40 % and particle size range 50-300 nm).
- ✓ PDW particle size is within number and intensity average DLS particle sizes (closer to d_p, N).
- ✓ Aggregation of polymer particles can be detected. Reduced scattering is more sensitive than particle size.
- ✓ The detection of new nucleations is challenging due to the big influence of the existing large polymer particles in the light scattering.

Acknowledgement



Thank you for
your attention



skal Herriko
ibertsitatea

IMAT



Universidad
del País Vasco
Euskal Herriko
Unibertsitatea

Inline monitoring of particle size in emulsion polymerization processes by Photon Density Wave (PDW) Spectroscopy

¹U.O. Aspiazu, ¹M. Paulis, ¹J.R. Leiza



NANOPAT:

¹ POLYMAT, Kimika Aplikatua saila, Kimika Fakultatea, University of the Basque Country UPV/EHU, Joxe Mari Korta zentroa, 20018 Donostia-San Sebastián.